Instruction manual

7-color LED Beauty Mask

Feature

- 1.100%Brand new and high quality
 - 2.Multiple functions:7 colors,manual adjustment and automatic transformation 3.Use it,you can stand,sit or lie down.
 - 4.Press and hold the start button for 3-5 seconds to turn on,and press and hold for 3-5 seconds to turn off.
 - 5.When starting the machine,click to switch to more colors.Automatically change color when you click the eighth time.
- Warm reminder:
- 1. Try to avoid looking directly at the light source of the mask, or looking at the light source for a long time. Use it with your eyes closed, or with a cotton pad or eye mask.
 - 2. Because the mask material is plastic, opening the package may have a smell, which is normal now, and the smell will soon disappear.
 - 3. The product can not be squeezed, can not be folded, cannot be heavy.
 - 4. Clean and luminous mask is strictly prohibited to soak in water, wipe it, elastic can be washed.
 - 5. Do not scrub the mask with corrosive chemical solvents.
 - 6. Avoid sun exposure.

Features:	Specification
Portable and Rechargeable	Item Type : LED Mask
One-buttonOperation	7 colors
Simple&Convenient	1*Main Machine
Long press:On/Off	1*Strap
Short press:Color Adjust	1*Manual
(7colors)	

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